

## I. Intended use

MagSi-DNA beads are ideal for purification or isolation of nucleic acids from various sources. The magnetic nanoparticles are intended as a solid phase extraction tool for custom buffer systems based on chaotropic as well non-chaotropic binding principles, and can be used for developing your own nucleic acid isolation and extraction methods, such as:

- Isolation of genomic, mitochondrial, or viral DNA from whole blood, cell lysates, human, animal, or plant tissue; isolation of RNA
- Isolation of genomic, plasmid, or phage DNA from bacterial cultures and bacteria from clinical samples (blood, stool, swabs, etc.)
- Clean-up of DNA from enzymatic reactions (restriction digestions, ligations) or chromatin immunoprecipitation (ChIP) procedures to remove excess primers, nucleotides, enzymes, salts, buffers and other substances that are unwanted in downstream applications

MagSi-DNA beads are magnetic silica beads with a highly dense magnetic core of iron oxide. Due to their magnetic properties, the beads typically collect within 10 seconds in a magnetic field. Their small particle size offers a very large active surface area and high binding capacity for nucleic acids. These particles are applicable in both manual and automated processes, *but continuous shaking* is needed because of their fast sedimentation.

## II. Principle

MagSi-DNA reversibly binds DNA and other nucleic acids under sample- and buffer-specific conditions. A solution containing DNA (e.g. lysate) is combined with the beads and an application-specific binding buffer. After incubation, nucleic acids are bound to the silica surface. By applying a suitable magnet to the container (tube/deepwell microplate) the bead pellet is separated from the sample mixture. Unwanted components are further removed by washing steps in a selection of buffers (alcohol/water solutions). Finally, nucleic acids are released in DNase/RNase-free water or buffer solution (e.g. Tris, Tris-EDTA, pH~8).

Silica and carboxylated (COOH) surfaces, but also nucleic acids, are negatively charged at neutral or basic pH, while both are also hydrated. For a chaotropic binding mechanism of DNA to particles, dehydration is needed. This can be achieved by for instance alcohol, and by agents such as guanidinium salts. Negative charges on the bead surface and the nucleic acid backbones are bridged by divalent cations. This can be reversed by a water solution.

For washing, mostly alcohol/water mixtures are used, which will keep the DNA in dehydrated form and bound to the beads. To reduce premature elution of DNA, salts can be added to the washing solution. Elution takes place in a low-salt conditions.

Non-chaotropic systems may use binding mechanisms with specific pH conditions, or binding by polyethylene glycol precipitation.

### Silica & Carboxylated MagSi-DNA beads

Optimal binding conditions differ for beads with silica or with carboxylated surfaces. In Table 1 below, some of the practical differences between the 2 types of beads are shown. (To develop a new application it is recommended to try both types in parallel! Contact magtivio for a test sample).

Table 1: Differences between silica (**MagSi-DNA**) and carboxylated (**MagSi-DNA COOH**) beads

Type of beads	Silica MagSi-DNA	Carboxylated MagSi-DNA COOH
<b>Compatible buffer systems</b>	Chaotropic buffers	PEG-based, low pH or chaotropic buffers
<b>Binding mechanism</b>	Precipitation with chaotropic salts	Precipitation by polymers like PEG, divalent cations (e.g. Mg <sup>2+</sup> ) or chaotropic salts
<b>Elution</b>	Low salt conditions	low salt conditions or pH shift from acidic binding to alkaline conditions

## III. Material supplied

- 2, 10, or 100 mL **MagSi-DNA** or **MagSi-DNA COOH** (supplied at 300 mg/mL, in storage buffer which is especially developed to keep the magnetic particles in a prolonged suspension time; containing 0.05% sodium azide )

### Additional materials needed

Depending on the application, reagents, equipment and consumables are needed:

- A specific set of lysis, binding, washing and elution buffers for the intended application
- Magnetic Separator for collection of the beads (see Order Information on the next page)
- Mixer/vortexer for homogenization of the beads and sample mixture.
- Optionally, a suspension buffer for preparation of the beads
- Container tubes or deep-well microplates and pipette tips

## IV. Product Usage

When stored at 2-8°C, this product is stable for up to 2 years, but no longer than the expiry date on the label. Store beads in well closed vial and in upright position to prevent drying of the beads since this makes them more difficult to re suspend. Do not freeze the product! Vortex bead suspension well before use.

MagSi-DNA (COOH) beads are suspended in sterile water with 0.05% sodium azide. The beads can be further pre-washed to avoid any impact in downstream applications. The suspension media can be replaced with your own buffer/storage media. The beads are compatible with typical organic solvents like ethanol or isopropanol. However, chemicals with strong redox-potential should be avoided.

The beads are stable in a pH range from 3 to 11 and at temperatures up to 95 degrees. After extensive incubations in these conditions, no degradation is detectable using spectrophotometric assays. Nevertheless, if you expect any interference in downstream applications, it is recommended to rinse the beads before use.


# MagSi-DNA

## MagSi-DNA COOH

### Product Description

Art.No.: MD0X017

Art.No.: MD0X019



Be aware that MagSi-DNA beads sediment within seconds and need homogenization during incubation. In some cases **MagSi-DNA 600, allround or 3.0** may be more suitable due to their longer sedimentation time, when separation speed is less critical.

#### IV. Protocols

The protocols below are intended as a guideline to develop a customized protocol and application.

##### A. Sample Preparation

Lyse your cell, tissue, or bacterial sample via:

- Using a surfactant like Tween 20/SDS/Triton X-100. Lysis efficiency may be improved by heating the sample mixture.
- mechanical disruption (sonication/French press)
- Enzymatic (lysozyme) methods

##### B. Binding

- Add the binding buffer of choice to the lysate and mix well to get a homogeneous suspension.
- Add beads. Mix beads by vortexing before adding them to the sample. Depending on the expected amount of DNA the volume of beads can be varied. A good starting point is 20  $\mu$ L when having 400-800  $\mu$ L of cell lysate.
- Mix sample and incubate 2-10 minutes to allow the DNA to bind to the bead surface.

##### C. Washing

- Following incubation, place sample tube in a magnetic separator.
- Wait until all the beads have been collected to the magnet. Discard the supernatant using a pipette, then remove the tube from the separator.
- Add wash buffer, vortex 10 seconds and place the sample tube in a magnetic separator in order to collect the beads and discard the supernatant.
- Wash the beads at least twice.

##### D. Elution

The Elution buffer consists of a nuclease-free, non-alcohol solution (TE-buffer) to rehydrate the DNA so it will elute from the bead. Concentrated TE-buffer can be added to the pure sample to improve storage properties.

- Elute DNA by adding 50-200  $\mu$ L elution buffer. Incubate 2-10 minutes at room temperature and mix several times.
- Collect beads with a magnetic separator and transfer the supernatant, containing the DNA, into a new tube.
- If eluate appears brown, repeat collection of the beads.
- Elution can be improved by repeating these steps or by incubating at 60°C during elution.

#### VI. Technical Data

Table 2: Technical data for MagSi-DNA and MagSi-DNA COOH

Product Name	MagSi-DNA (COOH)
Size	300 nm
Concentration	300 mg/mL
Magnetic content	95%
Surface Area	6.5 m <sup>2</sup> /g beads
Material	Magnetic silica beads optimized for nucleic acid isolation.
Solution additives	Water with 0.05% sodium azide
Storage	2-8°C

#### VII. Additional Information

##### Disclaimer

For Research Use Only (RUO). Not for drug, household or other uses. Safety Data Sheet (SDS) is available at [www.magtivio.com](http://www.magtivio.com).

##### Order Information

Product	Volume	Art. No.
MagSi-DNA	2 mL	MD01017
MagSi-DNA	10 mL	MD02017
MagSi-DNA	100 mL	MD03017

Product	Volume	Art. No.
MagSi-DNA COOH	2 mL	MD01019
MagSi-DNA COOH	10 mL	MD02019
MagSi-DNA COOH	100 mL	MD03019

##### Related Products

Product	Art. No.
MM-Separator M12 + 12	MD90001
MM-Separator M96	MD90002
MM-Separator 96 PCR	MDMG0005
MM-Separator 384 PCR	MDMG0006
MM-Separator 96 DeepWell	MDMG0013

##### magtivio B.V.

###### Office, Lab & Production

Daelderweg 9

6361 HK Nuth (The Netherlands)

Tel: +31-(0)45-208 4810

Fax: +31-(0)45-208 4817

E-mail: [info@magtivio.com](mailto:info@magtivio.com)